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Albumin ELISA Kit

Super sensitiv

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Albumin ELISA Kit

Super sensitive

For the in vitro determination of Albumin in urine and stool

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<u>Table of content</u>	<u>Page</u>
1. INTENDED USE	13
2. SUMMARY AND EXPLANATION OF THE TEST	13
3. PRINCIPLE OF THE TEST	13
4. MATERIAL SUPPLIED	14
5. MATERIAL REQUIRED BUT NOT SUPPLIED	14
6. PREPARATION AND STORAGE OF REAGENTS	15
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	16
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
TEST PROCEDURE	17
10. RESULTS	18
TYPICAL CALIBRATION CURVE	18
11. LIMITATIONS	19
12. QUALITY CONTROL	19
EXPECTED VALUES	19
13. PERFORMANCE CHARACTERISTICS	20
PRECISION AND REPRODUCIBILITY	20
SENSITIVITY	20
14. REFERENCES	20
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	21

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CE

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Albumin** in urine and stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Albumin reflects the major protein in human plasma (40-60%). It is synthesized in the liver depending of protein uptake. Albumin in fecal samples refers to an inflammation reaction combined with intestinal bleeding. Elevated levels of albumin and hemoglobin in stool are found not only by colorectal carcinomas, but also with polyps and during chronic inflammatory diseases (Morbus Crohn or Colitis Ulcerosa).

Indication

- Detection of source of bleeding in the lower gastrointestinal tract
- Detection of colorectal carcinoma
- Investigation of high risk patients
- Morbus Crohn, Colitis Ulcerosa

3. PRINCIPLE OF THE TEST

This Enzyme-Linked Immunosorbent Assay (ELISA) is a two step assay for the ultra sensitive determination of human Albumin in stool and urine. A polyclonal rabbit antibody specific for human Albumin is immobilized on microtiter plates and a second anti-Albumin antibody is conjugated to peroxidase (POD). Albumin in samples is bound to the immobilized antibody and after a washing step, to remove all unbound material the quantification of bound Albumin is carried out by adding a POD-labeled anti-Albumin antibody. After a second washing step substrate is added. The enzyme reaction is stopped and the absorbance is measured photometrically at 450 nm.

The intensity of the formed color is directly proportional to the concentration of Albumin in the samples.

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 6330MTP	precoated strips	12 x 8
K 6330WP	ELISA wash buffer concentrate (10x)	1 x 100 ml
K 6330PV	Sample dilution buffer, ready to use	1 x 100 ml
K 6330VP	Conjugate dilution buffer, ready-to-use	1 x 15 ml
K 6330K	Conjugate, (rabbit-anti-Albumin peroxidase-labeled)	1 x 50 µl
K 6330ST	Calibrators, lyophilized (0, 12.5, 50, 200, 800 µg/l)	4 x 5 vials
K 6330KO	Control, lyophilized	4 x 1 vial
K 6330TMB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 6330AC	ELISA stop solution, ready to use	1 x 7 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized water
- Precision pipettors calibrated to deliver 10-100 µl
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** should be diluted with aqua dest. **1:10** before use (add 450 ml aqua dest. to 50 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- **Calibrators** and **Controls** must be reconstituted with **500 µl aqua dest.** Allow the vial contents to dissolve for **10 minutes** and then mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted Calibrators and Controls are not stable and can not be stored.**
- The **Conjugate** (POD Antibody) must be diluted **1:400** in Conjugate dilution buffer (25 µl POD antibody in 10 ml buffer). The antibody is stable at 2-4°C until expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2-8°C.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Urine

Adjust urine to a pH 6 to 8 with 1 N NaOH. Samples can be stored for two weeks at 2-8 °C or at -20°C for longtime storage.

For the assay, dilute samples **1:200** with sample dilution buffer, e.g.

10 µl sample + 1990 µl sample dilution buffer.

Faeces

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the wash buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. 1 ml supernatant is given into an Eppendorf tube and centrifuged once more at 13.000 rpm for 5 min. For the assay, the supernatant is diluted **1:5** in sample dilution buffer (200 µl supernatant + 800 µl sample dilution buffer). **100 µl** of this solution is used for the assay.

We recommend to weight out the stool samples for each run new. Supernatant is not stable and can't be stored.

Stool samples can be stored at -20°C for 4 weeks. Avoid repeated thawing and freezing cycles.

Immundiagnostik recommends the use of the sample tubes from Roche Diagnostics/Mannheim (No. 745804) for sample preparation.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the precoated microtiter plate 5 x with 250 µl of diluted ELISA wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution. This is valid for all following washing steps. Perform each determination in duplicate for standards, controls and samples.


1. Add **100 µl** standard and control solutions and prediluted patient samples.
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
4. Add **100 µl** diluted Peroxidase-labeled anti-Albumin antibody.
5. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
7. Add **100 µl** TMB substrate solution. Incubate for **10-20 minutes** at room temperature.
8. Add **50 µl** stop solution and mix shortly.
9. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline-algorithm is recommended.

Typical calibration curve

Concentration [ng/ml]	800	200	50	12.5	0
OD mean values	1.353	1.114	0.687	0.269	0.118

The data is for demonstration only and cannot be used in place of data generations at the time of the assay.

Faeces

For the Albumin concentration of faeces samples, calculate as described in the following example:

Weight: 80 mg (1ml stool = 1g) = 0,08 ml

Dilution step 1: 5ml / 0,08ml = 62,5

Dilution step 2: 5

Dilution factor: 312.5

Multiply the result by **312.5** to get the real concentration. **The dilution factor depends on the weight of the faeces.**

Urine

To get the Albumin concentration in urine samples, the observed values must be multiplied by **200**.

11. LIMITATIONS

Samples with Albumin levels greater than the highest calibrator should be further diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

Normal range

Albumin in stool: < 9.2 mg/l (n = 76)

Albumin in urine: 5 - 16 mg/l

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Inter-Assay CV (n = 20)

Sample	Albumin [mg/l]	Inter-Assay CV [%]
1	50	8

Intra-Assay CV (n=20)

Sample	Albumin [mg/l]	Intra-Assay CV [%]
1	50	5

Sensitivity

The sensitivity or minimum detection limit of this assay is defined as the smallest single value that can be distinguished from zero at the 2 x standard variation as confidence limit. ($B0 + 2 SD = 12.5 \text{ ng/l}$)

Cross reactivity

No cross reactivity to other plasma proteins in faeces

14. REFERENCES

1. Träsch and Bloch: 1993; Clin. Lab. 39, 479
2. John et al.: 1994; Clin. Lab. 40, 77

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- For *in vitro* diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.