Chymotrypsin ELISA Kit

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For the in vitro determination of Chymotrypsin in stool

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Table of contents	Page
1. INTENDED USE	14
2. CLINICAL RELEVANCE	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
TEST PROCEDURE	18
10. RESULTS	20
TYPICAL CALIBRATION CURVE	20
11. LIMITATIONS	21
12. QUALITY CONTROL	21
EXPECTED VALUES	21
13. PERFORMANCE CHARACTERISTICS	22
SENSITIVITY	22
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Lagerung/Storage: 2 - 8 °C

1. INTENDED USE

The *Immundiagnostik* Assay is a sandwich Enzyme Immuno Assay intended for the quantitative determination of **Chymotrypsin** in stool. For in vitro diagnostic use only.

2. CLINICAL RELEVANCE

Chymotrypsin is a serine protease, which is secreted as excretory enzyme from the pancreas after food intake into the duodenum. Here food proteins are hydrolytically cleaved, preferentially next to aromatic residues. A small part of the active form of the enzyme is excreted within the stool. In the case of pancreas insufficiency secondary to a chronic pancreatitis, the secretion of the enzyme is reduced markedly, so that a reduced chymotrypsin activity can be measured in stool specimen.

The advantage of this assay is that the concentration of the enzyme is accurately determined and not the enzyme's activity.

Many studies have shown a longer stability of **chymotrypsin** (12 days at room temperature) compared to pancreatic elastase. Therefore, the **chymotrypsin ELISA** is a reliable and low-priced alternative to the pancreatic elastase assay.

Indications

- Chronic pancreatitis
- Exocrine pancreas insufficiency

3. PRINCIPLE OF THE TEST

This Enzyme Immuno Assay is a sandwich assay for the determination of Chymotrypsin in stool samples. The wells of the microtitre plate are coated with polyclonal antibodies directed against Chymotrypsin Protein. In a first incubation step, the Chymotrypsin in the samples is bound to the coated antibodies (in excess). To remove all unbound substances, a washing step is carried out.

In a second incubation step, a second antibody against Chymotrypsin is bound to the antigen. After another washing step, to remove all unbound substances, the solid phase is incubated with a peroxidase labelled Antibody. After a last washing step the plate is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stopping solution is then added. The

color converts to yellow. The intensity of the yellow color is directly proportional to the concentration of Chymotrypsin in the sample. A dose response curve of the absorbance (at 450 nm) unit vs. concentration is generated. Chymotrypsin, present in the patient samples, is determined directly from this calibration curve.

The combination of two specific antibodies in the Chymotrypsin ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

4. MATERIAL SUPPLIED

Catalogue No.	Kit Components	Quantity
K 6910MTP	one holder with precoated strips	96
K 6910WP	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6910A2	2 nd Antibody, (Mouse anti Chymotrypsin)	1 x 150 µl
K 6910K	Conjugate (POD Antikörper Goat anti Mouse)	1 x 150 µl
K 6910VP	2 nd Antibody dilution buffer	1 x 12 ml
K 6910ST	Calibrators, ready-to-use (1000; 250; 125; 62.5; 15.6; 0 ng/ml)	6 x 250 µl
K 6910ko	Control, ready-to-use	1 x 250 µl
K 6910TMB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 6910PV	Dilution buffer for the calibrators	1 x 11 ml
K 6910AC	ELISA stop solution, ready to use	1 x 7 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Distilled or deionized water
- Precision pipettes calibrated to deliver 10-1000 µl and disposable tips
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Absorbent paper
- Vortex-Mixer
- Microplate reader 450 nm

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are carefully stored at 2-8°C. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** should be diluted with aqua dest. **1:10** before use (add 450 ml aqua dest. to 50 ml ELISA wash buffer concentrate). Crystals may be formed due to high salt concentration. The crystals have to be dissolved **before dilution of the buffer concentrate** using a water bath (37°C). The buffer concentrates are stable at 2-8°C up to the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **2nd Antibody** has to be diluted **1:100** in 2nd antibody dilution buffer (100 µl 2nd antibody and 10 ml dilution buffer). The antibody is stable at 2 -8 °C up to the expiry date stated the label. **Diluted antibody solution is not stable and could not be stored.**
- The **POD labeled antibody** has to be diluted **1:100** in ELISA wash buffer (100 µl POD antibody and 10 ml ELISA wash buffer). The antibody is stable at 2 -8 °C up to the expiry date stated the label. **Diluted antibody solution is not stable and could not be stored.**

7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The stop solution consists of diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date stated on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Faeces

The test can be performed on either fresh or frozen stool samples. The samples should be refrigerated and could be stored at 2-8°C for 2 days. If the test cannot be performed within this period the specimen should be stored at -20°C or colder.

Add a stool sample of about 100 mg (size of a pea, please note the exact weight for the calculation) to 5 ml of the ELISA wash buffer and homogenize very thoroughly for 15 seconds on a Vortex-mixer. Centrifuge the suspension for 10 min at 3000 rpm. 1 ml of the supernatant is given into an Eppendorf tube and centrifuged once more at 13,000 rpm for 5 min. This supernatant could be stored at -20°C for about 1 month. 100 µl of this supernatant is used in the assay.

Immundiagnostik recommends the use of Roche Diagnostics / Mannheim sample preparation tubes, article No. 745804, for sample preparation.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

The precoated microtiter strips have to be washed **5 x with 250 µl** ELISA wash buffer before use. Bring all reagents to room temperature. Mix reagents well before use. Avoid direct sun light during all incubation steps. Covering the microtitre plate during the different incubation steps is recommended. Carrying out the tests in duplicate is also recommended.

1. Add **100 µl** calibrator dilution buffer into each well
2. Add 20 µl calibrators, control and samples and into each well.
3. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
4. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer.
5. Add **100 µl** prediluted 2nd antibody into each well.
6. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
7. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer.
8. Add **100 µl** prediluted Peroxidase-labeled antibody into each well.
9. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.


10. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
11. Add **100 µl** TMB substrate solution into each well.
12. Incubate for **5 - 15 minutes** at room temperature in the dark.
13. Add **50 µl** stop solution into each well and mix shortly.
14. Measure the extinction of the samples at **450 nm** directly after adding the stop solution.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

Typical calibration curve

Concentration [ng/ml]	1000	250	125	62.5	15.6	0
OD mean value	1.423	1.134	0.905	0.434	0.138	0.113

These data are for demonstration only and cannot be used instead of data obtained from the actual assay.

Faecal

To calculate the Chymotrypsin concentration of faecal specimen see the following example:

weight: 80 mg (1ml stool = 1g) = 0.08 ml

dilution step 1: 5 ml / 0,08 ml = 62.5

dilution factor: 62.5

Multiply the results with the calculated dilution factor (in this case 62.5) to get the Chymotrypsin concentration of the stool samples. **Please note:** the dilution factor depends on the weight of the used faecal specimen.

11. LIMITATIONS

Samples with Chymotrypsin levels greater than the highest calibrator, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

Stool: > 6.4 µg/ml mean made by Immundiagnostik (n=80)

These normal ranges should be used as a guideline only. It is recommended that each laboratory establishes an own expected range for its patient population.

13. PERFORMANCE CHARACTERISTICS

Sensitivity

n=20

Sample	Chymotrypsin mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.00625	0.00301	14.2

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for HVB and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any alterations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage.