EDN ELISA Kit

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
EDN ELISA Kit

For the in vitro Determination of EDN (Eosinophil- derived neurotoxin) in serum, plasma, urine and stool

valid from March 10th, 2005

K 6810

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of EDN (eosinophil-derived neurotoxin, eosinophil protein x, EPX) in serum, plasma, urine and stool. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

EDN (eosinophil-derived neurotoxin, eosinophil protein x, EPX), a cationic glycoprotein, which is released by activated eosinophiles, has strong cytotoxic characteristics and plays a significant role in the prevention of virus infections. It is released by the eosinophile granules in places where eosinophiles are mainly found, in the skin, lungs, urogenital, and gastrointestinal tract, that is, in the organs acting as an entry point for pathogens. The accumulation of EDN in the intestine is associated with inflammation and tissue damage.

Measuring of EDN in stool can serve as an objective parameter for a current clinical or sub-clinical chronic inflammation located in the gastrointestinal area. In the case of Colitis ulcerosa and Crohn's disease, the EDN measurement enables the evaluation of a disease's activity and the prediction of a relapse.

Indications

- Morbus Crohn
- Proof of a food allergy and incompatibility
- Assessment of an elimination diet
- Proof of damaged integrity of the intestinal mucous membrane caused by an invasive disease (e.g. CED, CC etc.)
- Proof of intestinal parasites
- Parasitoses

3 MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 6810MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6810WB	WASHBUF	ELISA wash buffer concentrate 10x	100 ml
K 6810PV	SAMPLEBUF	Sample dilution buffer 10x	100 ml
K 6810AP	ASYBUF	Assay buffer, ready-to-use	50 ml
K 6810ST	STD	Standard, ready-to-use (see specification for range)	7 vials
K 6810KO	CTRL	Control, ready-to-use (see specification for range)	1 vial
K6810K	CONJ	Conjugate, polyclonal peroxidase-labeled antibody, ready-to-use	12 ml
K 6810TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	20 ml
K 6810AC	STOP	ELISA stop solution, ready to use	20 ml

4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label. **Prepare just the appropriate amount necessary for the assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml a. bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be redissolved at room temperature or at 37°C using a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for one month.**
- The **sample dilution buffer** (SAMPLEBUF) should be diluted with aqua bidest. **1:10** before use (900 ml a. bidest. + 100 ml SAMPLEBUF), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be re-dissolved at room temperature or at 37°C using a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** could be stored in a closed flask at **2-8°C for one month.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

6. SAMPLE PREPARATION

Extraction step of the stool sample

The **sample dilution buffer** is used as sample extraction buffer. We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing of 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). Add 5 ml buffer to the stool sample in the preparation kit.

Constant buffer volume: 5 ml

Constant dilution factor: 1:50

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample independent of the sample amount.

Constant buffer volume: 5 ml

The dilution factor varies depending on the sample amount which has to be considered in the subsequent calculations as shown below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor is constant.

Variable buffer volume

Constant dilution factor: 1:50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results.

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix the weighed stool sample with the buffer and vortex for at least 30 sec. depending on the stool consistency.

Transfer approx. 1 ml stool suspension (**dilution step I**) in an Eppendorf-tube and centrifuge for 5 min at 13000 rpm (= 13000 g).

For analysis, pipette **100 µl** of the supernatant of **dilution step I** per well.

Dilution of samples

Urine samples

We recommend to analyze urine collected within 24 hours, whereby the EDN concentration is expressed as mg/day. If 24 h-urine sample is not available, urine from a single time point can be analyzed. In this case, the urinary creatinine should also be quantified, and the EDN results are presented as $\mu\text{g}/\text{mmol}$ creatinine.

Within 30 min of urine collection, the urine is separated by centrifugation, twice for 10 min at 1,350 x g and 4 °C. The supernatant is then transferred to a new plastic tube.

The urine samples should be diluted **1:51** with the assay buffer (ASYBUF), prior to analyses.

10 μl sample + 500 μl ASYBUF.

Serum/plasma samples

Fresh collected Serum/Plasma should be centrifuged within one hour. Store samples at -20 °C if not assayed on the same day. Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.

The **serum/plasma** samples should be diluted **1:5** with assay buffer (ASYBUF), prior to analyses.

50 μl sample + 200 μl ASYBUF.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes the two-site “sandwich” technique with two selected antibodies (monoclonal and polyclonal) that bind to human EDN.

Assay standards, controls and prediluted patient samples containing human EDN are added to wells of microplate that was coated with a high affine monoclonal anti-human EDN antibody. After the first incubation period, antibody immobilized on the wall of microtiter wells captures human EDN in the sample. Then a peroxidase-conjugated polyclonal anti-human EDN antibody is added to each microtiter well and a “sandwich” of capture antibody - human EDN – Peroxidase-conjugate is formed. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of EDN. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. EDN present in the patient samples, is determined directly from this curve.

Test procedure

<p>Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well</p>
<p>Mark the positions of STD /SAMPLE/CTRL (Standards/Sample/Control) on a protocol sheet</p>
<p>Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until expiry date stated on the label</p>
<p>Wash each well 5 times by dispensing 250 µl of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper</p>
<p>Add 100 µl of STD/SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective well. Use the corresponding sample dilution buffer as STD 0 ng/ml.</p>
<p>Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer</p>
<p>Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper</p>
<p>Add 100 µl CONJ (conjugate) into each well</p>

Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer

Discard the contents of each well. Wash 5 times by dispensing 250 µl of diluted WASHBUF (Wash buffer) into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper

Add 100 µl of SUB (substrate) into each well

Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark*

Add 100 µl of STOP (stop solution) into each well, mix thoroughly

Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Stool samples

To obtain the EDN concentration in stool samples, determined value by the dilution factor according to the sample preparation:

Sample preparation 1 and 2b: dilution factor constant : 1 : 50

Multiply the obtained result by **50** to get the final concentration.

Sample preparation 2a: dilution factor is variable.

The corresponding dilution factor is taken from the table and should be multiplied by the obtained result to get the final concentration.


Serum/plasma samples

For the calculation of the EDN concentration in **plasma/serum** the result has to be multiplied by **5**.

Urine samples

For the calculation of the EDN concentration in **urine** samples, the result has to be multiplied by 51.

Typical calibration curve

STD	OD1	OD2	mean OD	CV (%)
0	0.019	0.012	0.016	6.1
1	0.179	0.167	0.173	4.6
2	0.327	0.320	0.323	1.6
3	0.548	0.549	0.548	0.0
4	0.927	0.932	0.930	0.4
5	1.357	1.376	1.366	0.9
6	1.832	1.839	1.835	0.3
7	2.206	2.158	2.182	1.6

The data is for demonstration only and cannot be used for the evaluation of the user's results.

9. LIMITATIONS

Stool samples with EDN levels greater than the highest standard value, should be diluted with sample dilution buffer, and be re-assayed.

Serum/plasma and **urine** samples with EDN levels greater than the highest standard value, should be diluted with assay buffer, and be re-assayed.

10. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal ranges

(1 g stool is equivalent to 1ml)

Stool (n= 53): 357,6 ng/ml

Based on Immundiagnostik studies of evidently healthy persons (n= 53) a mean value of 357,64 ng/ml stool was estimated. For the present, the value of 1700 ng/ml stool should be considered as the upper limit of the test.

Urine (n = 50): 81,8 (26,7 – 164,2) µg/mmol Creatinine

Serum (n = 52): 26,4 (8,3 – 66,4) ng/ml

Plasma (n = 52): 18,1 (6,2 – 49,8) ng/ml

We recommend each laboratory to establish its own norm concentration range..

11. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

12. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG together with a written complaint.

06.10.2005 10032005_EDN.DOC