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Hb / Hp–Komplex ELISA Kit

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Table of contents	Page
1. INTENDED USE	18
2. SUMMARY AND EXPLANATION OF THE TEST	18
3. PRINCIPLE OF THE TEST	18
4. MATERIAL SUPPLIED	19
5. MATERIAL REQUIRED BUT NOT SUPPLIED	19
6. PREPARATION AND STORAGE OF REAGENTS	20
7. PRECAUTIONS	21
8. SPECIMEN COLLECTION AND PREPARATION	21
EXTRACTION OF THE STOOL SAMPLE	21
DILUTION OF SAMPLES	23
9. ASSAY PROCEDURE	23
PROCEDURAL NOTES	23
TEST PROCEDURE	24
10. RESULTS	25
TYPICAL CALIBRATION CURVE	25
11. LIMITATIONS	26
12. QUALITY CONTROL	26
EXPECTED VALUES	26
13. PERFORMANCE CHARACTERISTICS	27
PRECISION AND REPRODUCIBILITY	27
RECOVERY	28
SENSITIVITY	29
CROSS REACTIVITY	29
SAMPLE DILUTION	29
14. REFERENCES	29
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	30

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Manual

Hb / Hp-Complex ELISA Kit

*For the in vitro determination of Hemoglobin/Haptoglobin
Complex in stool*


Gültig ab/valid from 23.03.2006

K 7817

96

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Hemoglobin-/Haptoglobin-Complex** in stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

In contrast to other commercially available **Hemoglobin** quick tests, this **Hemoglobin-/Haptoglobin-Complex ELISA** does not require previous adherence to a diet (no raw meat etc.) and recognizes human **Hemoglobin** in 100-fold lower concentrations. This avoids false-negative results. Because of the choice of antibodies, false-positive results are almost excluded. Recent data show that the **Hemoglobin-/Haptoglobin-Complex** determination increases the clinical specificity and sensitivity.

Indications

- Occult blood in stool
- Crohn´s disease; Ulcerative Colitis
- Suspicion of colon carcinoma
- Polyps in the colorectum

3. PRINCIPLE OF THE TEST

The Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is used for quantitative determination of Hemoglobin-/Haptoglobin-Complex in stool.

Standards, controls and diluted patient samples which are assayed for Hemoglobin-/Haptoglobin-Complex are added to the wells of a microtiter plate coated with high affine anti-Haptoglobin antibodies. During the first incubation step, the Hemoglobin-/Haptoglobin-Complex is bound by the immobilized antibodies. A peroxidase-labeled antibody is then added into each microtiter well and a "sandwich" of capture antibody - Hemoglobin-/Haptoglobin-Complex – peroxidase-conjugate is formed. Tetramethylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the Hemoglobin-/Haptoglobin-Complex concentration of the sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standard. Hemoglobin-/Haptoglobin-Complex present in the patient samples is determined from this curve.

4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 7817MTP	PLATE	one holder with precoated strips	12 x 8
K 7817WP	WASHBUF	ELISA wash buffer concentrate (10x)	2 x 100 ml
K 7817PV	CONJ	Dilution buffer, ready-to-use	1 x 100 ml
K 7817K	SAMPLEBUF	Conjugate, (Rabbit anti Hb, POD-labeled)	1 x 200 µl
K 7817ST	STD	Standard, lyophil. (0; 0.1; 0.4; 1.6; 6.4 U/l)	5 vials
K 7817ko	CTRL	Controls	1 x 4 vials
K 7817TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 7817AC	STOP	ELISA stop solution, ready-to-use	1 x 7 ml

On request we will send you 3 calibrator sets free of charge. Any further calibrator sets will be charged.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution. The **buffer concentrate** (WASHBUF) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The lyophilized **standards** (STD) and **control** (CTRL) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards and control are **not stable**.
- The **conjugate** (CONJ) (POD-antibody) must be diluted **1:100** in **sample dilution buffer** (SAMPLEBUF) (100 µl CONJ + 10 ml SAMPLEBUF). The undiluted conjugate is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Extraction of the stool sample

Diluted wash buffer is used as a sample extraction buffer. We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample must be suspended in 5 ml wash buffer.

Constant buffer volume: 5 ml

Constant dilution factor: 1:50

2. Alternatively, stool samples can be manually weighted within the range of **80 – 120 mg**. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample not depending on the sample amount.

Constant buffer volume: 5 ml

The dilution factor varies depending on the sample amount which must be considered in the subsequent calculations as shown below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor remains constant.

Variable buffer volume

Constant dilution factor: 1:50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results.

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix stool sample and buffer; vortex for at least 30 sec. depending on the stool consistency.

Transfer approx. 1 ml stool suspension (**dilution step I**) to an Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm (= 13000 g). The resulting supernatant can be stored at **-20 °C for about 1 month**.

Dilution of samples

Stool samples

The supernatant of the centrifugation (dilution step I) is diluted **1:100 in sample dilution buffer (SAMPLEBUF)**. For example:

10 µl supernatant of dilution I + 990 µl SAMPLEBUF = 1:100 (**dilution step II**)

For analysis, pipette **100 µl** of the supernatant of **dilution step II** per well.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Stool samples can be stored at **-20°C** for 4 weeks. Avoid thawing and freezing cycles.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the precoated microtiter plate **5 x with 250 µl ELISA wash buffer**. Carry out the tests in duplicate. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution.


1. Add **100 µl STD/CTRL/SAMPLE** (Standards/Control/Sample) into respective well.
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** diluted WASHBUF.
4. Add **100 µl** pre-diluted **CONJ** (Peroxidase-labeled antibody).
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** diluted WASHBUF.
7. Add **100 µl SUB** (TMB substrate).
8. Incubate for **5-15 minutes** at room temperature.
9. Add **50 µl STOP** (stop solution) and mix shortly.
10. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL- algorithm is recommended.

Typical calibration curve

Concentration [U/I]	6.4	1.6	0.4	0.1	0
OD mean values	2.509	1.717	0.597	0.244	0.167

The data is for demonstration only and cannot be used for the evaluation of test results.

Faeces

Calculate Hemoglobin-/Haptoglobin-Complex concentration of stool samples as follows:

Sample weight: 80 mg (1ml stool = 1g) = 0,08 ml

Dilution step 1: 5 ml / 0,08 ml = 62,5

Dilution step 2: 100

Dilution factor: 62,5 x 100= 6250

For the calculation of the Hemoglobin-/Haptoglobin-Complex concentration the result must be multiplied by **6250**.

The dilution factor varies depending on the sample amount which must be considered in the subsequent calculation.

11. LIMITATIONS

Samples with Hemoglobin-/Haptoglobin-Complex levels greater than the highest calibrator value should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal range

Hemoglobin-/Haptoglobin-Complex (stool) < 2,5 U/ml (U/ml)

We recommend each laboratory to establish its own norm concentration range.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik Hemoglobin-/Haptoglobin-Complex (Hb/Hp) ELISA test was calculated from 20 replicate determinations on each of two samples.

Intra-Assay CV n= 20

Sample	Hb/Hp Mean value [mU/g]	Intra-Assay CV
		[%]
1	29.4	4
2	33.2	3

The total precision (inter-assay variation) of the Immundiagnostik Hemoglobin-/Haptoglobin-Complex (Hb/Hp) ELISA test was calculated from data on two samples obtained in 20 different assays by three technicians using two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

Sample	Hb/Hp Mean value [mU/g]	Inter-Assay CV
		[%]
1	39.2	8.7
2	38.7	9.7

Recovery

Two samples were spiked with Hemoglobin-/Haptoglobin calibrator and measured with this assay.

Recovery n=4

Sample [mU/g]	Spike [mU/g]	Hb/Hp Expected [mU/g]	Hb/Hp Measured
			[mU/g]
0.15	0.4	0.55	0.45
0.15	4	4.15	4.4
0.15	0.8	0.95	1.1
0.15	0.2	0.35	0.31
0.085	0.4	0.9	0.51
0.085	4	4.1	3.8
0.085	0.2	0.28	0.23
0.085	0.04	0.09	0.12

Sensitivity

The detection limit was defined as $B_0 + 2SD$. The zero-standard was measured 20 times.

n=20

Sample	Hb/Hp Mean value [OD]	Standard variation	Detection limit [U/ml]
1	0.060	0.005	0.15

Cross reactivity

No cross reactivity to other plasma proteins in stool.

Sample dilution

Linearity n= 2

Two patient samples were diluted with wash buffer. The results are shown below:

Sample	Dilution	Expected [mU/g]	Measured [mU/g]
A	undiluted	3.5	3.5
	1:2	1.7	1.5
	1:4	0.9	0.8
	1:8	0.4	0.4
B	undiluted	1.63	1.63
	1:2	0.8	0.5
	1:4	0.4	0.34
	1:8	0.2	0.2

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15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.