
# Helicobacter pylori Antigen ELISA Kit

*Zur in vitro Bestimmung des Helicobacter pylori in Stuhl*  
*For the in vitro Determination of Helicobacter pylori in stool*

Gültig ab / Valid from 22.03.2006

K 6920

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim  
Tel.: ++49 6251 70190-0  
Fax: ++ 49 6251 849430  
e.mail: [Info@immundiagnostik.com](mailto:Info@immundiagnostik.com)  
[www.Immundiagnostik.com](http://www.Immundiagnostik.com)

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# Helicobacter pylori Antigen ELISA Kit

*For the in vitro Determination of Helicobacter pylori in stool*


Valid from 22.03.2006

K 6920

96

## 1. INTENDED USE

This *Immundiagnostik Assay* is a sandwich ELISA for determination of *Helicobacter pylori* in stool. For research use only.

## 2. INTRODUCTION

*Helicobacter pylori* is regarded as a causative factor for chronic B-gastritis, drug-unrelated ulcer duodeni, and as an etiologic stimulus of gastric MALT-lymphoma. Furthermore, it is suspected of being involved in the pathogenesis of stomach carcinoma.

The epidemiology of infection by *Helicobacter pylori* has been characterized in western industrial nations by a linear increase with age increasing and in developing countries by a large number of children and juveniles being affected.

The currently used methods for the diagnosis of *Helicobacter pylori* infection are very sensitive and highly specific, but all include either invasive sampling or require special technical devices.

We have developed an ELISA for the detection of *Helicobacter pylori* infection from faecal samples. This cost-efficient methodology provides a reliable result without the loss of sensitivity or specificity and does not require invasive sampling.

### Indications:

- Detection of a *Helicobacter pylori* infection
- Monitoring the effect of a *Helicobacter pylori* treatment

## 3. PRINCIPLE OF THE TEST

A microtiter plate is coated with polyclonal antibodies specific for *Helicobacter pylori* which bind *Helicobacter pylori* from the patient sample in a first incubation step. After a washing step, the bound antigen, *Helicobacter pylori*, is incubated with a biotin-labelled antibody. After a further washing step, a peroxidase-labelled streptavidine is added. Then a chromogenic substrate (TMB) is added, followed by a stop solution. The colour changes from blue to yellow. The intensity of the colour is proportional to the amount of analyte (sample or control). The results are evaluated by comparison with a cut-off value.

## 4. MATERIAL SUPPLIED

Catalogue No.	Content	Kit Components	Quantity
K 6920MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6920WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6920PV	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 100 ml
K 6920VP	2.ABDIL	2. <sup>nd</sup> Antibody dilution buffer, ready-to-use	1 x 15 ml
K 6920KO1	CTRL NEG	Control negative, lyophilised	4 vials
K 6920KO2	CTRL POS	Control positive, lyophilised	4 vials
K 6920A2	2.AB	2. <sup>nd</sup> Antibody (anti-Helicobacter pylori biotin-labelled)	1 x 150 µl
K 6920K	CONJ	Conjugate, peroxidase-labelled streptavidine	1 x 150 µl
K 6920TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 6920AC	STOP	ELISA stop solution, ready to use	1 x 7 ml

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **2. antibody** (2.AB) must be diluted **1:100** in 2.<sup>nd</sup> antibody dilution buffer (2.ABDIL) (100 µl biotin-labelled antibody + 9.9 ml 2.<sup>nd</sup> antibody dilution buffer). The undiluted antibody is stable at 2 - 4 °C up to the expiry date stated on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **conjugate** (CONJ, POD-labelled streptavidine) must be diluted **1:100** in wash buffer (100 µl POD-streptavidine + 9.9 ml wash buffer). The undiluted conjugate is stable at 2-4 °C up to the expiry date stated on the label. **Diluted conjugate solution is not stable and can not be stored.**
- The **controls**, both **negative and positive** (CTRL NEG/CTRL POS), must be reconstituted with **500 µl** bidest. water. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted controls are not stable and can not be stored.**

## **7. SAMPLE PREPARATION**

### **Faeces**

The test can be performed on either fresh or frozen stool samples. If the test cannot be performed within one day, the specimen should be stored at  $-20^{\circ}\text{C}$  or colder.

Add a stool sample of **100 mg** to **1 ml** of the sample dilution buffer (SAMPLEBUF) and homogenize thoroughly on a Vortex-mixer. Centrifuge the suspension for 15 min at 3000 rpm.

*Immundiagnostik* recommends for sample preparation the use of Roche Diagnostics / Mannheim sample preparation tubes, article No. 745804.

## **8. ASSAY PROCEDURE**

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Carry out the tests in duplicate. Cover the microtiter plate during the different incubation steps. Bring all reagents samples to room temperature before use and mix well. Avoid direct sun light during all incubation steps. Wash the precoated microtiter strips **5 x with 250 µl** ELISA wash buffer before use. Afterwards dry the microtiter plate on absorbent paper.

1. Add **100 µl** of **CTRL POS/CTRL NEG/ SAMPLEBUF/SAMPLE** (control positive, control negative, sample dilution buffer (blank) and samples (stool suspension)) in duplicate into respective well.
2. Cover and incubate for **1 hour**, shaking on a horizontal mixer at room temperature.
3. Aspirate and wash the wells **7 x with 250 µl** diluted WASHBUF (Wash buffer).
4. Add **100 µl** diluted **2.AB** (2<sup>nd</sup> antibody, biotin-labelled) into each well.
5. Cover and incubate for **1 hour**, shaking on a horizontal mixer at room temperature.
6. Aspirate and wash the wells **7 x with 250 µl** diluted WASHBUF (Wash buffer).
7. Add **100 µl** diluted **CONJ** (Conjugate, Peroxidase-labelled) into each well.
8. Cover and incubate for **1 hour**, shaking on a horizontal mixer at room temperature.
9. Aspirate and wash the wells **7 x with 250 µl** diluted WASHBUF (Washbuffer).
10. Add **100 µl SUB** (substrate) into each well.
11. Incubate for **5-15 minutes** at room temperature in the dark.
12. Add **50 µl STOP** (stop solution) into each well and mix shortly.
13. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the positive control exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## 9. RESULTS

### *Calculation of the values*

Blank absorbance must be  $< 0.150$ . If the absorbance of the blank is higher, the analysis must be repeated.

### *Test result*

**Cut off-value** = 0,150 at 450/620 nm and 0,190 at 450 nm

Samples with absorbance more than 0,020 (measurement at 450/620 nm) or 0,025 (measurement at 450 nm) above the cut off-value are considered as **positive**.

Samples with an absorbance 0,020 (0,025) below or above the cut off-value are regarded as **border line** samples and must be re-analyzed.

If the value from the repeated measurement is again in the **border line** range, a new sample should be analyzed.

If the value from the repeated measurement is again in the range above the cut off-value, the result is considered as **positive**.

Samples with absorbance more than 0,020 (0,025) under the cut off-value are classified as **negative**.

## 10. QUALITY CONTROL

The test is correct when the absorption of the negative control at 450 nm is less than 0.190 and at 450/620 nm less than 0.150, and the positive control value at 450 nm greater than 0.54 and at 450/620 nm greater than 0.50. If the observed values differ from the expected, the test must be repeated.

Immundiagnostik AG recommends the use of commercial available controls for internal quality control.

**Control samples** should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods.

## **11. PRECAUTIONS**

- For research use only.
- Quality control guidelines should be observed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

## **12. TECHNICAL HINTS**

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiry date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

### **13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

24/03/2006 HELICO\_30.06.2005.DOC