Hämoglobin ELISA Kit

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Table of contents	Page
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	17
Extraction of the stool sample	17
Dilution of samples	19
9. ASSAY PROCEDURE	19
Procedural notes	19
Test procedure	20
10. RESULTS	21
Typical calibration curve	21
11. LIMITATIONS	22
12. QUALITY CONTROL	22
Expected values	22
13. PERFORMANCE CHARACTERISTICS	22
Sensitivity	22
14. REFERENCES	23
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Manual

Hemoglobin ELISA Kit

For the in vitro determination of Hemoglobin in stool


Gültig ab/valid from 31.03.2006

K 7816

96

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of human **Hemoglobin** in stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

In contrast to other commercially available **Hemoglobin** quick tests, this **Hemoglobin ELISA** does not require previous adherence to a diet (no raw meat etc.) and recognizes human **Hemoglobin** in 100-fold lower concentrations. This avoids false-negative results. Because of the choice of antibodies, false-positive results are almost excluded.

Recent data show, that the clinical specificity and sensitivity can be increased by using **Hemoglobin ELISA**.

Indication

- Occult blood in stool
- Crohn´s disease; Ulcerative Colitis
- Suspicion of colon carcinoma
- Polyps in the colorectum

3. PRINCIPLE OF THE TEST

The Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is used for quantitative determination of Hemoglobin in stool. The Hemoglobin in the sample is bound to anti-Hemoglobin antibodies (in excess), which are immobilized on the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step an anti-Hemoglobin peroxidase labeled antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Hemoglobin in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the results obtained from the calibrators. Hemoglobin, present in the patient samples, is determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 7816MTP	PLATE	one holder with pre-coated strips	12 x 8
K 7816WP	WASHBUF	ELISA wash buffer concentrate (10x)	2 x 100 ml
K 7816PV	SAMPLEBUF	Dilution buffer, ready-to-use	1 x 100 ml
K 7816K	CONJ	Conjugate, (sheep-anti Hb, peroxidase-labeled)	1 x 200 µl
K 7816ST	STD	Calibrators, lyophil. (100; 33; 11; 3.3; 0 ng/ml)	5 vials
K 7816KO	CTRL	Control, lyophil.	1 vial
K 7816TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 7816AC	STOP	ELISA stop solution, ready-to-use	1 x 7 ml

On request, we will send you 3 calibrator sets free of charge. Any further calibrator sets will be charged.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution. The **buffer concentrate** (WASHBUF) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The lyophilized **standards** (STD) and **control** (CTRL) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards and control are **not stable**.
- The **conjugate** (CONJ) (POD-antibody) must be diluted **1:100** in **sample dilution buffer** (SAMPLEBUF) (100 µl CONJ + 10 ml SAMPLEBUF). The undiluted conjugate is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Extraction of the stool sample

Diluted wash buffer is used as a sample extraction buffer. We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample must be suspended in 5 ml wash buffer.

Constant buffer volume: 5 ml

Constant dilution factor: 1:50

2. Alternatively, stool samples can be manually weighted within the range of **80 – 120 mg**. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample not depending on the sample amount.

Constant buffer volume: 5 ml

The dilution factor varies depending on the sample amount which must be considered in the subsequent calculations as shown below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor remains constant.

Variable buffer volume

Constant dilution factor: 1:50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results.

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix stool sample and buffer; vortex for at least 30 sec. depending on the stool consistency.

Transfer approx. 1 ml stool suspension (**dilution step I**) to an Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm (= 13000 g). The resulting supernatant can be stored at -20 °C for about 1 month.

Dilution of samples

Stool samples

The supernatant of the centrifugation (dilution step I) is diluted **1:10 in sample dilution buffer (SAMPLEBUF)**. For example:

100 µl supernatant of dilution I + 900 µl SAMPLEBUF = 1:10 (**dilution step II**)

For analysis, pipette **100 µl** of the supernatant of **dilution step II** per well.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Stool samples can be stored at -20°C for 4 weeks. Avoid thawing and freezing cycles.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any variations of the test procedure, without consulting the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.


1. Add **100 µl STD/CTRL/SAMPLE** (Standards/Control/Sample) into respective well.
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** diluted WASHBUF.
4. Add **100 µl** pre-diluted **CONJ** (Peroxidase-labeled antibody).
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** diluted WASHBUF.
7. Add **100 µl SUB** (TMB substrate).
8. Incubate for **10-20 minutes** at room temperature.
9. Add **50 µl STOP** (stop solution) and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

10. RESULTS

A calibration curve is constructed from the standards. For data evaluation, commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL- algorithm is recommended.

Typical calibration curve

Concentration [ng/ml]	100	33	11	3.3	0
OD mean value	1.667	0.792	0.334	0.156	0.073

The data is for demonstration only and cannot be used for the evaluation of test results.

Faeces

For the Hemoglobin concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 10

dilution factor: 62,5 x 10= 625

Multiply the result with **625** to get the real concentration. The dilution factor is dependent on the weight of the faeces.

11. LIMITATIONS

Samples with Hemoglobin levels greater than the highest calibrator, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

Normal ranges

Hemoglobin (stool)	<10 µg/ml
	Grey area 5-10 µg/ml

13. PERFORMANCE CHARACTERISTICS

Sensitivity

The sensitivity was set as $B_0 + 2SD$. The zero-standard was measured 20 times. A detection limit of 1.8 ng/ml was obtained.

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15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.