# Lysozym ELISA Kit

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# Lysozyme ELISA Kit

*For the in vitro determination of Lysozyme in serum and stool*

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Table of content	Page
<b>1. INTENDED USE</b>	<b>17</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>17</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>17</b>
<b>4. MATERIAL SUPPLIED</b>	<b>18</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>19</b>
<b>7. PRECAUTIONS</b>	<b>20</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>20</b>
SERUM	20
FAECES	20
<b>9. ASSAY PROCEDURE</b>	<b>21</b>
PROCEDURAL NOTES	21
TEST PROCEDURE	21
<b>10. RESULTS</b>	<b>23</b>
TYPICAL CALIBRATION CURVE	23
<b>11. LIMITATIONS</b>	<b>24</b>
<b>12. QUALITY CONTROL</b>	<b>24</b>
EXPECTED VALUES	24
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>25</b>
PRECISION AND REPRODUCIBILITY	25
RECOVERY	26
SENSITIVITY	26
MEAN VALUE	26
STANDARD	26
DETECTION LIMIT	26
CROSS REACTIVITY	27
SAMPLE DILUTION	27
<b>14. REFERENCES</b>	<b>27</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>28</b>

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Lagerung/Storage: 2 - 8 °C

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Lysozyme** in serum and stool. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

**Lysozyme** (muramidase) is a protein with a molecular weight of approx. 15 kD and belongs to the group of alkaline glycosidases. **Lysozyme** is made by granulocytes, monocytes and macrophages. The main source for faecal **Lysozyme** are the intestinal granulocytes. **Lysozyme** can be detected in all cells of the inflammatory infiltrate during an acute attack of Crohn's disease. To some extent, **Lysozyme** is also secreted actively by mononuclear cells into the bowel lumen.

### Indication

- Diagnosis and monitoring of Crohn's Disease, Boech's Disease (in serum)
- Bacterial, viral, allergenic or autoimmune caused bowel inflammations

## 3. PRINCIPLE OF THE TEST

This Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is suitable for the quantitative determination of Lysozyme from human serum and stool. In a first incubation step, the Lysozyme in the patient samples is captured by the precoated antibody. The quantification of the bound Lysozyme is carried out by adding an rabbit anti Lysozym- Peroxidase-labelled conjugate. Finally the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is used to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Lysozyme in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. Lysozyme, present in the patient samples, is determined directly from this curve.

#### 4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 6900MTP	one holder with precoated strips	96
K 6900WB	ELISA wash concentrate 10x	2 x 100 ml
K 6900VP	POD dilution buffer	1 x 15 ml
K 6900K	Conjugate, (Sheep-anti-Lysozym, Peroxidase-labelled)	1 x 50 $\mu$ l
K 6750ST	Calibrators, ready to use (0; 1.1; 3.3; 10; 30 ng/ml)	5 x 1 ml
K 6750KO	Control	1 x 1 ml
K 6750TMB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 6750AC	ELISA stop solution, ready to use	1 x 7 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized water
- Precision pipettors calibrated to deliver 50-100  $\mu$ l
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **washing buffer concentrate** should be diluted with Aqua dest. **1:10** before use (add 450 ml Aqua. dest. to 50 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **Conjugate** (POD antibody) has to be diluted **1:1000** in POD dilution buffer (10 µl POD antibody and 10 ml POD dilution buffer). The antibody is stable at 2-8 °C until expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2-8°C.

## 7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop Solution consists of diluted Sulfuric Acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### *Serum*

Samples are diluted between **1:500** and **1:1000**.

Use this dilution factor to calculate the Lysozyme concentration.

### *Faeces*

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the washing buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. 1 ml supernatant is given into an eppendorf tube and centrifuged once more at 13.000 rpm for 5 min. The resulting supernatant could be stored at -20°C for about 1 month.

After thawing, the supernatant is centrifuged at 13.000 rpm for 2 min. Afterwards the supernatant is diluted **1:2** in washing buffer (200µl supernatant is added to 200 µl sample dilution buffer). **100 µl** of this solution is used for the assay.

Immundiagnostik recommends the use of the sample tubes from Boehringer / Mannheim for sample preparation.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** calibrators and patient samples (serum and stool diluted, see above).
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the cavities **5 x with 250 µl** of the diluted ELISA wash buffer solution.
4. Add **100 µl** of Peroxidase-labelled conjugate.
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the cavities **5 x with 250 µl** of the diluted washing buffer solution.
7. Add **100 µl** of TMB substrate solution.
8. Incubate for **10-20 minutes** at room temperature until colour differences are sufficient.


9. Add **50 µl** of stop solution and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## 10. RESULTS

A calibration curve is constructed from the calibrators. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

Concentration [ng/ml]	30	10	3.3	1.1	0
OD mean value	1.257	1.05	0.519	0.224	0.063

The data is for demonstration only and cannot be used in place of data generations at the time of the assay.

### Serum

The result has to be multiplied with **500** or **1000** to calculate the serum value.

### Faeces

For the Lysozyme concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 2

dilution factor:  $62,5 \times 2 = 125$

Multiply the result with **125** to get the real concentration. The dilution factor depends on the weight of the faeces.

## 11. LIMITATIONS

Samples with Lysozyme levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### *Expected values*

Lysozyme concentration (faeces): < 600 ng/ml

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

The precision (intra-assay variation) of the Immundiagnostik Lysozyme ELISA test was calculated from 20 replicate determinations on each of one samples.

Intra-Assay CV n= 20

Sample	Lysozyme [ng/ml]	Intra-Assay CV [%]
1	6.8	8
2	1.4	13

The total precision (inter-assay variation) of the Immundiagnostik Lysozyme ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

Sample	Lysozyme [ng/ml]	Inter-Assay CV [%]
1	5.6	17

*Recovery*

Two samples were spiked with Lysozyme calibrator and measured with this assay.

Recovery n=2

Sample [ng/ml]	Spike [ng/ml]	Lysozyme expected [ng/ml]	Lysozyme measured [ng/ml]
1.5	1	2.5	3.5
1.5	3	4.5	5.5
1.5	6	7.5	7.1
0.75	15	15.8	12.1
0.75	5	5.8	5.7
0.75	1.7	2.4	2.7
0.75	0.7	1.4	1.5

*Sensitivity*

n=20

Sample	Lysozyme mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.184	0.029	0.5

### *Cross reactivity*

No cross reactivity to other proteins in serum and stool.

### *Sample dilution*

Linearity n= 2

Two patient serum samples were diluted with ELISA wash buffer. The results are shown below:

Sample	Dilution	Expected [ng/ml]	Measured [ng/ml]
A	undiluted	15.3	15.3
	1:2	7.7	7.5
	1:4	3.8	3.9
	1:8	1.9	<detection limit
B	undiluted	9.3	9.3
	1:2	4.65	6.0
	1:4	2.3	4.1
	1:8	1.16	1.6

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## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from this.