# MPO ELISA Kit

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*For the in vitro determination of Myeloperoxidase (MPO) in stool  
and urine*

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<b>1. INTENDED USE</b>	<b>17</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>17</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>17</b>
<b>4. MATERIAL SUPPLIED</b>	<b>18</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>19</b>
<b>7. PRECAUTIONS</b>	<b>20</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>20</b>
<b>9. ASSAY PROCEDURE</b>	<b>21</b>
PROCEDURAL NOTES	21
TEST PROCEDURE	21
<b>10. RESULTS</b>	<b>23</b>
TYPICAL CALIBRATION CURVE	23
<b>11. LIMITATIONS</b>	<b>24</b>
<b>12. QUALITY CONTROL</b>	<b>24</b>
EXPECTED VALUES	24
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>25</b>
PRECISION AND REPRODUCIBILITY	25
RECOVERY	26
SENSITIVITY	26
SAMPLE DILUTION	27
CROSS REACTIVITY	27
<b>14. REFERENCES</b>	<b>27</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>28</b>

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•  
*For the in vitro determination of Myeloperoxidase (MPO) in stool  
and urine*

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Lagerung/Storage: 2 - 8 °C

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Myeloperoxidase** in urine and stool, designed to be also suitable for small series of specimen. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

The granules of neutrophils (approx. 70% of the white blood cells) contain a large number of different enzymes. **Myeloperoxidase** (MPO) catalyzes the oxidation of substances through  $H_2O_2$ . The **MPO**  $H_2O_2$ -system has a toxic effect on many micro-organisms such as bacteria, fungi, viruses and mycoplasma. The efficiency of the bacteria-destructive Myeloperoxidase  $H_2O_2$ -system is increased by PMN-Elastase. **MPO** determination in the stool reflects the inflammatory activity of Crohn's disease or ulcerative colitis.

### Indication

- Marker for inflammatory activities in the gastrointestinal tract
- Renal transplant rejection
- Oxidative stress
- For the differentiation between allergic and infectious asthma

## 3. PRINCIPLE OF THE TEST

This Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is suitable for the quantitative determination of Myeloperoxidase in urine and stool. In a first incubation step, the Myeloperoxidase in the samples is bound to an available excess of antibodies against Myeloperoxidase, which are immobilised to the surface of the microtitre plates. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labelled antibody against MPO is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of MPO in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. MPO, present in the patient samples, is determined directly from this curve.

#### 4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 6630MTP	one holder with precoated strips	96
K 6630WP	ELISA wash buffer concentrate 10x	1 x 100 ml
K 6630ST	Standarts, lyophilized (see specification for range)	4 x 5 vials
K 6630KO	Control, lyophilized (see specification for range)	4 x 1 vial
K 6630A2	Detection antibody, (mouse monoclonal anti-MPO antibody, biotinylized), concentrate	1 x 200 µl
K 6630K	Conjugate, (streptavidin peroxidase labeled), concentrate	1 x 200 µl
K 6630PV	Sample dilution buffer, ready to use	1 x 50 ml
K 6630TMB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 6630AC	ELISA stop solution, ready to use	1 x 7 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water
- Precision pipettes calibrated to deliver 10-100 µl
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Horizontal mixer
- Vortex
- Microplate reader 450 nm

## 6. PREPARATION AND STORAGE OF REAGENTS

- The **ELISA wash buffer concentrate** should be diluted with aqua dest. **1:10** before use (add 900 ml a. dest. to 100 ml concentrate). Crystals may be formed due to high salt concentration. The crystals have to be dissolved **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C up to the expiry date stated on the label. Diluted solutions can be stored at 2-8°C for 1 month.
- The **standards and controls** have to be reconstituted with **500 µl** aqua dest. Allow the vial to stand for 10 minutes and then mix thoroughly by gentle conversion to insure complete reconstitution. Reconstituted standards and controls are **not stable**.
- The **detection antibody (AB)** has to be diluted **1:101** in wash buffer (WASHBUF) (100 µl AB + 10 ml WASHBUF). The antibody is stable at 2-8 °C until expiry date given on the label. **Diluted antibody solution is not stable and could not be stored.**
- The **conjugate (CONJ)** has to be diluted **1:101** in wash buffer (WASHBUF) (100 µl CONJ and 10 ml WASHBUF). The antibody is stable at 2-8 °C until expiry date given on the label. **Diluted conjugate is not stable and could not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2-8°C. The antibody coated microtitre plate should be stored dry and dark at 2-8°C.

## 7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The stop Solution consists diluted sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date stated on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Faeces

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the wash buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. 1 ml supernatant is given into an Eppendorf tube and centrifuged once more at 13.000 rpm for 5 min. The resulting supernatant could be stored at -20°C for about 1 month.

After thawing, the supernatant is centrifuged at 13.000 rpm for 2 min. Afterwards the supernatant is diluted **1:10** in wash buffer (100 µl supernatant is added to 900 µl wash buffer). **100 µl** of this solution is used for the assay.

Immundiagnostik recommends the use of the sample tubes from Boehringer / Mannheim for sample preparation.

### Urine

Urine have to stored at -20°C. For assay use dilute the urine samples **1:10** in sample dilution buffer (100µl urine in 900µl sample dilution buffer).

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- Carry out the assay with the actual manual delivered with the kit.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** standard, control and prediluted patient samples.
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of biotin-labelled detection antibody solution.
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
7. Add **100 µl** of Peroxidase-labelled antibody conjugate solution.
8. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
9. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.


10. Add **100 µl** of substrate solution (TMB substrate solution).
11. Incubate for **5-15** minutes at room temperature until colour differences are sufficient.
12. Add 50 µl of STOP (stop solution) into each well, mix thoroughly.
13. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

## 10. RESULTS

A calibration curve is constructed from the calibrators measured in the assay. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

Concentration [ng/ml]	100	33	11	3.6	0
OD mean value	1.322	0.543	0.232	0.106	0.048

These data are for demonstration only and cannot be used instead of data obtained from the actual assay.

## Faeces

For the MPO concentration of faeces samples, calculate as described in the following example:

Amount of sample: 80 mg (1ml stool = 1g) = 0,08 ml

Dilution step 1: 5ml / 0,08ml = 62,5

Dilution step 2: 10

Dilutionfactor: 62,5 x 10 = 625

The from the calibration curve calculated MPO concentration have to be multiplied by **625** to get the right concentration. **Caution: the multiplication factor is dependend on the amount of stool, one use at the beginning.**

## Urine samples:

The from the calibration curve calculated MPO concentration have to be multiplied by **10** to get the right concentration.

## 11. LIMITATIONS

Samples with Myeloperoxidase levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

**Immundiagnostik recommends commercial control samples for internal quality control.**

Control samples or pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### *Expected values*

Myeloperoxidase concentration

Stool: < 2000 ng/g

These normal ranges should be used as a guideline only. It is recommended that each laboratory establishes an own expected range for its patient population.

### 13. PERFORMANCE CHARACTERISTICS

#### *Precision and reproducibility*

The precision (intra-assay variation) of the Immundiagnostik MPO ELISA test was calculated from 20 replicate determinations on each of one samples.

Intra-Assay CV n= 20

Sample	MPO mean value [ng/ml]	Intra-Assay CV [%]
1	147.1	4.3
2	288.6	4.8

The total precision (inter-assay variation) of the Immundiagnostik MPO ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

Sample	MPO mean value [ng/ml]	Inter-Assay CV [%]
1	171.7	12
2	239.9	15

*Recovery*

Two samples were spiked with MPO and measured with this assay.

Recovery n=2

Sample [ng/ml]	Spike [ng/ml]	MPO expected [ng/ml]	MPO measured [ng/ml]
116	500	616	514
116	320	436	401
116	200	316	336
116	125	241	254
92	500	592	504
92	320	412	388
92	200	292	297
92	125	217	204

*Sensitivity*

n=20

Sample	MPO mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.013	0.003	1.6

### Sample dilution

Linearity n= 2

Two patient serum samples were diluted with wash buffer. The results are shown below:

Sample	Dilution	Expected [ng/ml]	Measured [ng/ml]
A	1:40	14.5	14.5
	1:80	7.2	7.1
	1:160	3.6	3.5
B	1:40	19.5	19.5
	1:80	9.75	10.1
	1:160	4.8	5.2

### Cross reactivity

No cross reactivity to other plasma proteins in stool.

Alpha-1-Antitrypsin 0 %

Albumin 0 %

CRP 0 %

Lysozyme 0 %

slgA 0 %

PMN-Elastase 0 %

Calprotectin 0 %

## 14. REFERENCES

1. Saiki T et al.: 1998; Kurume Med. J. 45, 69

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for HBV and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any alterations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable responsible for any damage.