Vitamin A/E

HPLC

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For the determination of Vitamin A/E in plasma and serum

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1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	16
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
SAMPLE AND STANDARD PREPARATION	17
CHROMATOGRAPHIC CONDITIONS	18
10. TREATMENT OF THE COLUMN	18
11. RESULTS	18
CALCULATION	18
TYPICAL CHROMATOGRAM	19
12. LIMITATIONS	19
13. QUALITY CONTROL	19
EXPECTED VALUES	19
CONTROLS	19
14. PERFORMANCE CHARACTERISTICS	20
PRECISION AND REPRODUCIBILITY	20
LINEARITY	20
DETECTION LIMIT	20
RECOVERY	20
15. DISPOSAL	21
16. TROUBLESHOOTING	21
17. REFERENCES	22
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Comstock G.W., Alberg A.J., Helzlsouer K.J. (1993). Reported effects of long-term freezer storage on concentrations of retinol, β -carotene, and α -tocopherol in serum or plasma summarized. Clin Chem 39/6; 1075-1078

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Vitamin A/E

HPLC – Application

*For the determination of Vitamin A/E in
plasma and serum*

Catalogue No.: KC 1600
Package size: 100 determinations
Storage: 20-25 °C standard and int. standard -20 °C

CE

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of Vitamin A/E in serum and plasma. This Assay is designed for in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin A (Retinol) and Vitamin E (Tocopherol) are fatsoluble vitamins, which could be stored for longer periods in the adipose tissue. Both, lack and excess could express in complaints.

Vitamin A is essential for the visual process and recovers the skin and mucosa. A lack of Vitamin A will reduce the visual power up to total blindness. An excess of Vitamin A could cause headache, damage of the skin, liver disease, painful alteration in the skeleton or foetal damage.

Vitamin E (Tocopherol) is an antioxidant and protects unsaturated fatty acids against oxidation. It also protects the cells of the body by catching radicals.

A lack of vitamin E in animal experiments demonstrates diseases of muscle, nervous system, heart, liver and reproduction system. These symptoms were not observed in humans. Vitamin E can be stored in the adipose tissue in large amounts. A lack can be caused by a malfunction in digestion or resorption of fatty acids.

3. PRINCIPLE OF THE TEST

The first step in the determination of Vitamin A and E includes the sample preparation. In the first step an internal standard solution is added. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is used for injection into the HPLC system.

The separation via HPLC follows an isocratic method at 25°C using a „reversed phase“ column; one run lasts 15 minutes. The detection is performed by an UV detector at two different wavelengths (Vitamin A: 325 nm; Vitamin E: 300 nm). The quantification is performed with the delivered standard solution; the concentration is calculated via integration of the peak areas in the internal standard calibration mode.

Summary:

The application of Vitamin A and E for HPLC allows the determination of both vitamins in an easy, fast, and precise method. The kit includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As many other parameters the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test series the one-point calibration is much more economic than 6-point calibration for immuno assays).

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
KC1600lm	Mobile phase	1000 ml
KC1600st	Standard (conc. is given on the label)	10 ml
KC1600is	Internal Standard	5 ml
KC1600fr	Precipitation reagent	50 ml
KC1600vl	Dilution solution	10 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with UV-detector
- Reversed phase C₁₈-column

6. PREPARATION AND STORAGE OF REAGENTS

All reagents are ready to use and become delivered in soluble form.

All reagents are stable at 20-25 °C, the standard solution and the internal standard at -20 °C up to the date of expiry (see label of the test package).

7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum and EDTA-plasma are suited as samples.

The sample is light- and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately.

The samples are stable in the dark at 2-8°C for 12 h (Vitamin A) and 48 h (Vitamin E). At -20°C Vitamin A is stable for 1 month, Vitamin E for 3 month.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

Sample and standard preparation

Pipette into 1.5 ml reaction tubes:

Standard:

250 µl standard solution

+

50 µl internal standard solution

+

250 µl dilution solution

+

250 µl precipitation solution

Samples and controls:

250 µl sample or control

+

50 µl internal standard solution

+

500 µl precipitation solution

Mix well. Leave the tubes for **30 minutes at 2-8°C** and centrifuge afterwards at 10.000g for 10 minutes.

Inject **100 µl** of the supernatant for chromatography into the HPLC-system

Chromatographic conditions

Column material:	Nucleosil C ₁₈ ; 10 µm; 150 x 4 mm
Flow rate:	0.8 ml/min
Detection:	UV
	Vitamin A: 325 nm
	Vitamin E: 300 nm

The wavelength should be switched after 7 min

Running time per chromatogram: 15 minutes

To avoid contamination of the next run, mobile phase has to be used as autosampler wash solution.

Immundiagnostik recommends to use a guard-column to enlarge lifetime of the column.

10. TREATMENT OF THE COLUMN

After analysis the column could be left in the mobile phase. Before use, the system should be equilibrated with ca. 50 ml eluent.


11. RESULTS

Calculation

$$\frac{\text{peak area patient} \cdot \text{concentration of the standard}}{\text{peak area Internal standard patient}} * F = \text{concentration patient sample}$$

$$F = \frac{\text{peak area Internal standard of the standard - mixture}}{\text{peak area of each vitamin of the standard - mixture}}$$

Typical chromatogram

12. LIMITATIONS

Hemolytic and lipemic samples should not be measured

13. QUALITY CONTROL

Expected values

Vitamin A: 200 - 800 µg/l

Vitamin E: 3 - 14 mg/l

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate

statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay VK:	Vitamin A:	1,9 % (0,55 mg/l)	[n=6]
	Vitamin A:	1,2 % (1,18 mg/l)	[n=6]
	Vitamin E:	1,5 % (9,0 mg/l)	[n=6]
	Vitamin E:	1,1 % (14,9 mg/l)	[n=6]
Inter-Assay VK:	Vitamin A:	4,9 % (0,6 mg/l)	[n = 8]
	Vitamin A:	3,1 % (1,0 mg/l)	[n = 8]
	Vitamin E:	4,6 % (8,3 mg/l)	[n = 8]
	Vitamin E:	4,7 % (24 mg/l)	[n = 8]

Linearity

Vitamin A:	up to 10 mg/l
Vitamin E:	up to 50 mg/l

Detection limit

Vitamin A:	0.05 mg/l
Vitamin E:	1 mg/l

Recovery

Vitamin A:	98.8 %
Vitamin E:	101 %

15. DISPOSAL

The mobile phase, precipitation reagent, internal standard and standard must be disposed as non-halogenated solvent.

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signalcord and connection
	Detectorlamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Doublepeaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min

Problem	Possible reason	Solution
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not calm	Pump delivers imprecise	Check pump, degas the system
	Detector flowcell is dirty	Clean flow cell

17. REFERENCES

Sushil K.J., Mc Coy B., Wise R. (1994). Vitamin E and the hypercoagulability of neonatal blood. Clin Cim Acta 225; 97-103.

Comstock G.W., Alberg A.J., Helzlsouer K.J. (1993). Reported effects of long-term freezer storage on concentrations of retinol, β -carotene, and α -tocopherol in serum or plasma summarized. Clin Chem 39/6; 1075-1078.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.