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Vitamin B₁ HPLC Kit

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*For the determination of Vitamin B₁(thiamin pyrophosphate) in
EDTA-blood*

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1. INTENDED USE	15
2. SUMMARY AND EXPLANATION OF THE TEST	15
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	17
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	17
7. PRECAUTIONS	18
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
SAMPLE AND STANDARD PREPARATION	19
CHROMATOGRAPHIC CONDITIONS	20
10. TREATMENT OF THE COLUMN	20
11. RESULTS	20
CALCULATION	20
TYPICAL CHROMATOGRAM	21
12. LIMITATIONS	21
13. QUALITY CONTROL	21
EXPECTED VALUES	21
CONTROLS	22
14. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
LINEARITY	22
DETECTION LIMIT	22
15. DISPOSAL	22
16. TROUBLESHOOTING	23
17. REFERENCES	24
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24

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Lagerung/Storage: 4 °C, Calibrator and controls at -20 °C

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of Vitamin B1 in EDTA-blood. This assay is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin B₁ (thiamin) is a watersoluble Vitamin, which consists of a pyrimidine- and thiazolring linked via a methylen bridge. It is sensitive against alkaline solution, oxidation and reduction.

Vitamin B₁ is produced by plants and microorganisms. It is found free, peptide bound and as phospho-esters (mono-, di- and triphospho-esters). In animals and also in humans it is necessary to be supplemented by food. The intake of thiamin in the gut is maintained by active transport and passive diffusion. The different phospho-esters are synthesized by phosphorylation and dephosphorylation. The active form in metabolism is thiamin pyrophosphate and thiamin triphosphate in brain. After dephosphorylation thiamin is secreted by the kidney.

Thiamin pyrophosphate plays an important role as a co-enzyme in carbohydrate- and amino acid metabolism. An important reaction is the oxidative carboxylation. Thiamin itself is required for stimulating nerve cells. Beside this, it stimulates the fatty acid- and cholesterol-synthesis in nervous tissues.

A classical disease for the lack of Vitamin B₁ is Beri Beri, which is known from asians eating predominantly white rice. The symptoms are paralysis, drop in muscle mass and heart failure. Other diseases are the Wernicke-encephalopathie, the Korsakow-syndrome and several forms of the Landry's paralysis. Also many alcoholics have deficient thiamin status.

Applications:

- determination of Vitamin B₁ status
- disturbance in amino acid metabolism
- malabsorption (alcoholism)
- suspicion for neuritis

3. PRINCIPLE OF THE TEST

The first step in the determination of thiamin pyrophosphate includes the sample preparation with additional derivatisation. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is used for the derivatisation (10 min at 60°C), Vitamin B1 will be transformed into a fluorescent product. The sample is cooled, centrifuged and injected into the HPLC system.

The separation via HPLC follows an isocratic method at 25°C using a „reversed phase“ column; one run lasts 15 minutes. The quantification is performed with the calibrator delivered; the concentration is calculated via integration of the peak areas.

Summary

The application of vitamin B1 for HPLC allows the determination of thiamin pyrophosphate in an easy, fast, and precise method. The kits includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As with many other parameters the advantage of HPLC measurements are the simultaneous handling of many analytes in a single test. The „complete HPLC-system“ enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test runs the one-point calibration is much more economic than 6-point calibration for immuno assays).

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
KC 2201LM	Mobile phase	1 x 1000 ml
KC 2201KA	Calibrator, lyophilized	1 vial
KC 2201FR	Precipitating reagent	5 ml
KC 2201VL	Dilution solution	20 ml
KC 2201RB	Reaction buffer	5 ml
KC 2201DL	Derivatisation solution (lyoph.)	5.5 ml
KC 2201LC	Solution C	5.5 ml
KC 2201KO	Control 1 and 2; 250 µl lyophilized	2 x 3 vials

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C₁₈-column
- Heatable shaker or waterbath
- Vortex mixer

6. PREPARATION AND STORAGE OF REAGENTS

- The **calibrator**, EDTA- whole blood with a defined concentration of thiamin pyrophosphate should be resuspended with dilution solution (volume is given on the label), aliquoted and stored at -20 °C. The concentration of Vitamin B1 is varying slightly from batch to batch; the actual concentration is given on the label.

- Reconstitute the **controls** in **250 µl** dilution solution.
- Resuspend the **derivatisation solution** in 5.5 ml solution C. The dissolved derivatisation solution is stable for 6 month at 2-8 °C.
- All other test components are ready to use.
- Store test reagents at 2-8 °C, the **calibrator and the controls** at -20 °C. They are stable up to the expiry date which is given on the label.

7. PRECAUTIONS

- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The supplied reagents contain reagents like acetonitrile (mobile phase) and acid (precipitating reagent). Although diluted, it still must be handled with care and should only be handled with gloves, eye protection, and appropriate protective clothing. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA-whole-blood has to be used in this assay.

Vitamin B1 is light- and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately.

The samples are stable in the dark at 2-8°C for 1 day. For longer storage samples should be frozen at -20°C. Do not freeze the samples a second time.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.

- The assay should always be performed due to the manual which is given in the kit.

Sample and standard preparation

Pipette into 1.5 ml reaction tubes:

50 µl sample (EDTA-whole blood), calibrator or control

+

150 µl dilution solution

+

50 µl precipitating reagent

Mix well. Leave the tubes for **10 minutes at 2-8°C** and centrifuge afterwards at 10.000 x g for 10 minutes.

Take **150 µl** supernatant

+

50 µl reaction buffer

+

50 µl derivatisation solution,

Incubate for **10 minutes at 60°C** on a shaker or in a waterbath.

Cool to 2-8°C. Centrifuge at 10.000 x g for 5 minutes.

Inject **50 µl** of the supernatant for chromatography into the HPLC-system

The prepared supernatant is stable for 3 days at 2-8 °C.

Chromatographic conditions

Column material:	Bischoff Prontosil Eurobond; 5 µm; 125 x 4 mm MERCK Lichrospher RP 18; 125 x 4 mm
Flow rate:	0.8 - 1.2 ml/min
Detection:	Fluorescence excitation: 365 nm emission: 440 nm
Running time per chromatogram:	12 minutes

10. TREATMENT OF THE COLUMN


After the analysis the column should be flushed with 50 ml aqua bidest (1 ml/min) and stored in 20% methanol in aqua bidest (ca. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 50 ml mobile phase.

11. RESULTS

Calculation

$$c_{B_1 \text{ patient}} = \frac{\text{Peak height patient} \times \text{conc. calibrator (ng/ml)}}{\text{Peak height calibrator}}$$

Typical chromatogram

12. LIMITATIONS

We recommend not to measure lipaemic patient samples. The measurement of serum and plasma samples is possible but not recommended, because the concentration is mostly below the detection limit.

13. QUALITY CONTROL

Expected values

Vitamin B1 normal values: 32 – 95 ng/ml (mean \pm 2SD)

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or EDTA-blood pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay VK: 3,3 % (31,2 ng/ml) [n = 6]
4,3 % (59,0 ng/ml) [n=6]

Inter-Assay VK: 3,2 % (33,0 ng/ml) [n = 12]
4,7 % (62,3 ng/ml) [n = 12]

Linearity

up to 250 ng/ml

Detection limit

0.5 ng/ml

15. DISPOSAL

The mobile phase and derivatisation solution must be disposed as non-halogenated solvent. The precipitation solution could be neutralized with NaOH and if the pH value is neutral it can be disposed as salt solution. (**Important:** Reaction will produce heat, be careful)

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

17. REFERENCES

1. Tallaksen C.M.E., T. Bohmer, H. Bell (1991). Concomitant determination of thiamin and its phosphate esters in human blood and serum by high-performance liquid chromatography. *J. Chromatogr*, 564, 127-136.
2. Herbeth B., J. Zittoun, L. Miravet, M. Boureay-Causse, G. Carre-Guery, E. Delacoux, C. Le Devehat, A. Lemoine, J.P. Mareschi, J. Martin, G. Portier de Courcy and J. Sancho (1986). Reference intervals for vitamin B1, B2, E, D, retinol, and folate in blood: Usefulness of dietary selection criteria. *Clin. Chem.* 32/9, 1756-1759.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.