Vitamin B₆ HPLC Kit

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For the determination of Vitamin B₆ in plasma and serum

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| | |
|---|-----------|
| 1. INTENDED USE | 14 |
| 2. SUMMARY AND EXPLANATION OF THE TEST | 14 |
| 3. PRINCIPLE OF THE TEST | 14 |
| 4. MATERIAL SUPPLIED | 15 |
| 5. MATERIAL REQUIRED BUT NOT SUPPLIED | 16 |
| 6. PREPARATION AND STORAGE OF REAGENTS | 16 |
| 7. PRECAUTIONS | 16 |
| 8. SPECIMEN COLLECTION AND PREPARATION | 17 |
| 9. ASSAY PROCEDURE | 17 |
| PROCEDURAL NOTES | 17 |
| SAMPLE AND STANDARD PREPARATION | 18 |
| CHROMATOGRAPHIC CONDITIONS | 18 |
| 10. TREATMENT OF THE COLUMN | 19 |
| 11. RESULTS | 19 |
| CALCULATION | 19 |
| TYPICAL CHROMATOGRAM | 19 |
| 12. LIMITATIONS | 20 |
| 13. QUALITY CONTROL | 20 |
| EXPECTED VALUES | 20 |
| CONTROLS | 20 |
| 14. PERFORMANCE CHARACTERISTICS | 21 |
| PRECISION AND REPRODUCIBILITY | 21 |
| LINEARITY | 21 |
| DETECTION LIMIT | 21 |
| RECOVERY | 21 |
| 15. DISPOSAL | 21 |
| 16. TROUBLESHOOTING | 22 |
| 17. REFERENCES | 23 |
| 18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE | 23 |

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1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of vitamin B6 in plasma and serum. This Assay is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin B6 is the term for pyridoxin, pyridoxal and pyridoxamin and the appropriate phosphate products. All forms can be transformed into the active form pyridoxal-5-phosphate (here described as „vitamin B6“).

Vitamin B6 functions as a coenzyme and is essential for more than 50 reactions in the protein metabolism thereby synthesizing, transforming or degrading amino acids. Vitamin B6 supports the resorption of amino acids and their transport into the cells. Furtheron vitamin B6 contributes to the synthesis of neuro transmitters and amine products (histamin).

Due to the fact that vitamin B6 contributes to a variety of different reactions lack of vitamin B6 results in various clinical pictures.

- Applications:
- determination of vitamin B6 status
 - homocysteinaemia
 - skin diseases
 - movement disorders
 - Anaemia, depression

3. PRINCIPLE OF THE TEST

The first step in the determination of Vitamin B6 includes the sample preparation with additional derivatisation. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is used for derivatisation (20 min at 60°C) thereby transforming the vitamin B6 into a fluorescent product. The sample is cooled, centrifuged and injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a „reversed phase“ column; one run lasts 10 minutes. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak areas.

Summary:

The application of vitamin B6 for HPLC makes it possible to determine the vitamin in an easy, fast and precise method. The kits includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As many other parameters the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test series the one-point calibration is much more economic than 6-point calibration for immuno assays).

4. MATERIAL SUPPLIED

| Catalogue No | Kit Components | Quantity |
|--------------|-------------------------------------|-------------|
| KC 2100LM | Mobile phase | 1 x 1000 ml |
| KC 2100KA | Calibrator, 4 ml lyophilized | 1 vial |
| KC 2100DL | Derivatisation solution | 3 x 8.5 ml |
| KC 2100FR | Precipitating reagent | 1 x 5 ml |
| KC 2100RE | Reconstitution solution | 1 x 10 ml |
| KC 2100KO | Control 1 and 2, 250 µl lyophilized | 2 x 3 vials |

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C₁₈-column - Bischoff ProntoSil Eurobond C18 5 µm, 125 x 4 mm
- Heatable shaker or water bath
- Vortex mixer

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** in **4 ml** reconstitution solution. Take aliquots and store them at -20°C. Avoid repeated freeze-thaw circles. The concentration of vitamin B6 might have minor changes from lot to lot.
- Reconstitute the **controls** in **250 µl** reconstitution solution.
- All other test reagents are stable at 2-8 °C, up to the date of expiry stated on the label.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The precipitating reagent contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.

- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA-plasma and serum is suited for this test system.

Vitamin B6 is light- and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately.

The samples are stable in the dark at 2-8°C for 1 week. For longer storage samples should be frozen at -20°C.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

10. TREATMENT OF THE COLUMN


After the analysis the column should be flushed with 50 ml aqua bidest (1 ml/min) and stored in 20% methanol in aqua bidest (ca. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 50 ml Eluent.

11. RESULTS

Calculation

$$c_{B_6 \text{ patient}} = \frac{\text{Peak height patient} \times \text{conc. calibrator (ng / mL)}}{\text{Peak height calibrator}}$$

Typical chromatogram

12. LIMITATIONS

To minor interferences whole blood should be diluted 1:1 with aqua bidest. The calculated concentration must then be multiplied with 2.

13. QUALITY CONTROL

Expected values

Vitamin B6 normal values: 4.3 – 17.9 ng/ml (n=90)

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV: 2.8 % (24.7 ng/ml) [n = 12]

Inter-Assay CV: 4.1 % (24.6 ng/ml) [n = 12]

Linearity

up to 250 ng/ml

Detection limit

0.2 ng/ml

Recovery

95,6 %

15. DISPOSAL

The derivatisation solution can be oxidized with hydrogen peroxide and after the pH value is adjusted inbetween 6-8, it can be disposed as aqueous salt solution. The mobile phase and the precipitation solution could be neutralized with NaOH and if the pH value is neutral it can be disposed as salt solution. (**Important:** Reaction will produce heat, be careful)

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

| Problem | Possible reason | Solution |
|--------------------------|---|---|
| No signal | No or defect connection to evaluation system | Check signalcord and connection |
| | Detectorlamp is altered | Change lamp |
| No peaks | Injector is congested | Check Injector |
| Doublepeaks | Dead volume in fittings and / or column | Renew fittings and / or column |
| Contaminating peaks | Injector dirty | Clean injector |
| | Contamination at the head of the column | Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase |
| | Air in the system | Degas pump |
| | Autosampler vials contaminated | Use new vials or clean them with methanol |
| Broad peaks, tailing | Precolumn / column exhausted | Use new precolumn / column |
| Variable retention times | Drift in temperature | Use a column oven |
| | Pump delivers imprecise | Check pump, degas the system |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |
| Baseline is drifting | Detector lamp did not reach working temperature yet | Wait |
| | Detector lamp is too old | Renew lamp |
| | System is not in steady state yet | Rinse system mobile phase for 15 min |
| | Pump delivers imprecise | Check pump, degas the system |
| Baseline is not smooth | Pump delivers imprecise | Check pump, degas the system |
| | Detector flowcell is dirty | Clean flow cell |

17. REFERENCES

Ambrosch A. et al. (2000). Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic Med* 18; 185-192.

Dierkes J. et al. (2001). Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis* 158; 161-164.

Dierkes et al. (2001). Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *J Renal Nut* 11; 67-72.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.