Vitamin C HPLC Kit

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For the determination of Vitamin C in Li-heparine plasma

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1. INTENDED USE	14
2. SUMMARY AND EXPLANATION	14
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
SAMPLE AND STANDARD PREPARATION	17
CHROMATOGRAPHIC CONDITIONS	18
10. TREATMENT OF THE COLUMN	18
11. RESULTS	18
CALCULATION	18
TYPICAL CHROMATOGRAM	19
12. LIMITATIONS	19
13. QUALITY CONTROL	19
EXPECTED VALUES	19
CONTROLS	20
14. PERFORMANCE CHARACTERISTICS	20
PRECISION AND REPRODUCIBILITY	20
LINEARITY	20
DETECTION LIMIT	20
15. DISPOSAL	20
16. TROUBLESHOOTING	21
17. REFERENCES	22
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	22

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1. INTENDED USE

The Immundiagnostik HPLC application is a chromatographic procedure intended for the quantitative determination of Vitamin C in plasma and serum. This application is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION

In the 15th to 17th century more sailors died of scorbout than of any other disease. The food provided on board contained nearly no Vitamin C. In the 16th century the importance of Vitamin C, supplied by citrus fruits in healing scorbout was discovered.

Ascorbic acid (Vitamin C) is a strong reducing substance. The oxidation of Vitamin C leads over a radical intermediate to dehydroascorbic acid, *in vivo*. The three forms mentioned constitute a reversible redox-system.

Ascorbic acid plays an important role in hydroxylation reactions, i.e. in the synthesis of collagen. So it is rather important for the *de novo* synthesis of bone, cartilage and tooth and for the healing of wounds. Vitamin C is needed for the production of noradrenalin. Another important role of Vitamin C is that of an antioxidant, that means it protects other substances from the oxidizing effects of oxygen. Ascorbic acid promotes the resorption of iron in the intestine. It reduces the production of nitrosamines which might cause cancer.

The primary unspecific signals of a lack of Vitamin C are: tiredness, physical and mental weakness and increased susceptibility for infections. Psychic disturbances like depressions or hysteria are possible.

The advantage of HPLC analysis is the simultaneous handling of many analytes in a single test. The complete HPLC system enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routine diagnostics. In most cases a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per assay. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled.

Indication

- Determination of Vitamin C status

3. PRINCIPLE OF THE TEST

The application of Vitamin C for HPLC allows the determination of Vitamin C in an easy, fast, and precise method. The kit includes all reagents in ready to use for preparation and separation of the samples. The columns are not included.

The first step in the determination of Vitamin C is a fast and easy sample preparation. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is injected into the HPLC system.

The separation via HPLC follows an isocratic method at 25°C using a „reversed phase“ column; one run lasts 12 minutes. Vitamin C is detected by an UV-monitor. The quantification is performed with the delivered calibrator; the concentration is calculated via integration of the peak areas by the external standard method.

4. MATERIAL SUPPLIED

Artikel Nr.	Kit Komponenten	Menge
KC2900LM	Mobile phase	1000 ml
KC2900KA	calibrator (0,25 ml; lyoph. concentration is given on the label)	8 vials
KC2900FR	Precipitating reagent (lyoph.)	1 vial
KC2900RL	Reconstitution solution	27 ml
KC2900KO	Control 1 and 2; 250 µl lyophilized (concentration, see the product data sheet)	2 x 3 vials

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Vortex- mixer
- Centrifugation tubes (1.5 ml)
- High precision pipettes
- HPLC pump with UV-detector
- reversed phase C₁₈-column
- Centrifuge

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** with 250 µl aqua bidest. The reconstituted standard solution is not stable and can not be used for a second time. Do **not freeze** it again. The concentration is given on the label.
- Reconstitute the **controls** in 250 µl aqua bidest.
- Reconstitute the **precipitating reagent** in 25 ml reconstitution solution in a ultrasonic bath (approx. 10 min.). The precipitation solution is stable for 3 months at 2-8 °C.
- All reagents are stable at 2-8 °C, **calibrator and controls** at -20 °C up to the date of expiry (see label of the test package).

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The precipitating reagent contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Venous fasting blood is suitable for this test system. We recommend to use Lithium-heparine plasma because of a better stability of vitamin C. Therefore commercial available sample tubes (e.g. Sarstedt S-Monovette LH) should be used. The best stability is achieved when a 7.5 ml tube is filled with 2 ml of blood.

Vitamin C is highly sensitive against oxidation; therefore samples should be stabilized immediately after arrival in the laboratory. For stabilisation, the precipitating reagent is added (see 9. Assay procedure).

Plasma or serum, containing the precipitating reagent is stable for 24 h at 2-8°C. After centrifugation the supernatant is stable for 3 month at -20°C.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

Sample and standard preparation

Add into 1.5 ml reaction tubes:

200 µl patient sample, standard and controls

+

200 µl precipitating reagent

Mix well. Incubate the tubes for **10 minutes at 2-8°C** and then centrifuge at 10.000 x g for 10 min.

Inject **20 µl** of the supernatant (The supernatant is stable at least for 24 hours at room temperature, if kept in the dark).

Chromatographic conditions

Column material:	Bischoff Prontosil AQ 5 µm
Dimension:	125 x 4 mm
Flow rate:	0.75 ml/min
UV-Detector:	254 nm
Injection volume:	20 µl
Running time/sample:	12 min

10. TREATMENT OF THE COLUMN


After each run the column should be washed with 50 ml aqua bidest (1 ml/min) and stored in 20% methanol (v/v in aqua bidest, 30 ml), flow rate 0.7 ml/min. Before use, the system should be equilibrated with ca. 50 ml eluent.

11. RESULTS

Calculation

$$C_{patient} = \frac{\text{Peak area patient} \times \text{conc. calibrator (mg/l)}}{\text{Peak area calibrator}}$$

Typical chromatogram

12. LIMITATIONS

EDTA-blood is not suitable for this test system and should not be used.

13. QUALITY CONTROL

Expected values

4 to 20 mg/l (mean \pm 2 SD).

It is recommended, that each laboratory should develop its own normal range. The values mentioned above are established in Immundiagnostik AG laboratories and could deviate from other published data.

Controls

Control samples or plasma pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intraassay CV%: 1.6 (10.3 mg/l) [n=10]

Interassay CV%: 2.7 (7.3 mg/l) [n=10]

Linearity

up to 250 mg/l

Detection limit

0.25 mg/l

15. DISPOSAL

The mobile phase and the precipitation solution could be neutralized with NaOH and if the pH value is neutral it can be disposed as salt solution. (**Important:** Reaction will produce heat, be careful)

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

17. REFERENCES

1. Hultqvist M. et al. (1997). Plasma concentrations of vitamin C, vitamin E and/or malondialdehyde as markers of oxygen free radical production during hemodialysis. Clin Nephrol 47; 37-46.
2. Falch J.A. (1998). Low levels of serum ascorbic acid in elderly patients with hip fracture. Scand J Clin Lab Invest 58; 225-228.
3. Ballmer et al. (1994). Depletion of plasma vitamin C but not of vitamin E in response to cardiac operations. J Thorac Cardiovasc Surg 108; 311-320.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.