## ELISA Kit

# Anti Gewebetransglutaminase IgA

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18.01.2006 tTG IgA\_10012005.DOC

# Anti Tissue Transglutaminase IgA

*For the in vitro determination of anti-tissue-Transglutaminase  
IgA antibodies (anti-tTG) in serum and plasma*

Valid from 06.09.2005

Catalogue No.: K 9399

Package size: 96 Determinations

Storage: 2 - 8 °C

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of human **anti-tissue-Transglutaminase-IgA antibodies** in serum and plasma. For *in vitro* diagnostic use only.

## 2. PRINCIPLE OF THE TEST

This Enzyme Immuno Assay is a sandwich assay for the determination of anti-tissue transglutaminase antibodies (anti-tTG IgA) in serum samples. The wells of the microtiter plate are coated with the antigen. In a first incubation step, the anti-tTG antibodies are bound to the immobilized antigen. To remove all unbound substances, a washing step is carried out. Then a peroxidase-labeled antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the peroxidase substrate, Tetramethylbenzidine (TMB). An acidic stopping solution is then added. The color converts to yellow. The results are evaluated using a cut off value.

## 3. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 9399MTP	One holder with precoated strips	12 x 8 wells
K 9399WP	ELISA wash buffer concentrate 10x	2 x 100 ml
K 9399K01	Control negative, lyophilized	4 vials
K 9399K02	Control positive, lyophilized	4 vials
K 9399K	POD antibodies (rabbit-anti-IgA antibody, peroxidase labeled)	1 x 200 µl
K 9399TMB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 9399AC	ELISA stop solution, ready-to-use	1 x 7 ml

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 5-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100  $\mu$ l** should be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **POD antibody** must be diluted **1:100** in ELISA wash buffer (e.g. 100  $\mu$ l POD antibody + 10 ml ELISA wash buffer). The antibody is stable at 2-8 °C until expiry date stated on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **Controls** (positive and negative) must be reconstituted with **500  $\mu$ l washing buffer**. Allow the vial to stand for 10 minutes and then mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted **controls are not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 6. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 7. SPECIMEN COLLECTION AND PREPARATION

### Serum/Plasma

Dilute serum and plasma **1:250** in ELISA wash buffer in two dilution steps.

Example:

**Dilution step 1: 10 µl serum + 90 µl wash buffer**

**Dilution step 2: 40 µl dilution 1 + 960 µl wash buffer**

## 8. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Performe assay always according the enclosed manual.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** control solutions and prediluted patient samples per well.
2. Incubate for **1 hour** at room temperature (21-26°C).
3. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** of prediluted, peroxidase-labeled anti-IgA antibody into each well.
5. Incubate for **1 hour** at room temperature, shaking on a horizontal mixer.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
7. Add **100 µl** of TMB-substrate solution into each well.
8. Incubate for **10-20 minutes** at room temperature.
9. Add **50 µl** of stop solution into each well and mix shortly.
10. Determine absorption immediately with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

## 9. RESULTS

### Calculation

$$\begin{aligned} & \text{mean OD negative reference sera} * \text{Factor(See data sheet)} \\ & = \text{cut off} = 18 \text{ AU/ml} \end{aligned}$$

### Example for Cut off:

$$0.104 \times 2.5 \text{ (factor data sheet)} = 0.260 \text{ (cut off)} = 18 \text{ AU/ml}$$

The calculated cut off OD is 0.260.

Samples with a mean OD higher than the negative reference sera are considered as positive:

### Example for positive patient sample

Patient sample extinction = 0.685

$$0.520 = 18 \text{ AU/ml}$$

$$0.685 = X \qquad X = \frac{0.685 * 18}{0.260} = 47.4 \text{ AU/ml}$$

## 10. LIMITATIONS

Samples with h tTG IgA levels greater than the positive control value should be diluted and re-assayed.

## 11. QUALITY CONTROL

Immundiagnostik recommends commercially control samples as internal quality control.

**Control samples** should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

htTG IgA Serum cut off > 22 AU/ml positiv

htTG IgA Serum cut off < 16 AU/ml negativ

Results between 16 - 22 AU/ml are in the grey range.

We recommend each laboratory to establish its own normal range.

## 12. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

1/18/2006 tTG IgA10012004.DOC