# slgA ELISA Kit

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# slgA ELISA Kit

*For the in vitro determination of Secretory IgA in saliva and stool*

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Table of content	Page
<b>1. INTENDED USE</b>	<b>16</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>16</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>17</b>
<b>4. MATERIAL SUPPLIED</b>	<b>17</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>18</b>
<b>7. PRECAUTIONS</b>	<b>19</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>19</b>
<b>9. ASSAY PROCEDURE</b>	<b>20</b>
PROCEDURAL NOTES	20
TEST PROCEDURE	20
<b>10. RESULTS</b>	<b>22</b>
TYPICAL CALIBRATION CURVE	22
<b>11. LIMITATIONS</b>	<b>23</b>
<b>12. QUALITY CONTROL</b>	<b>23</b>
EXPECTED VALUES	23
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>24</b>
PRECISION AND REPRODUCIBILITY	24
RECOVERY	25
SENSITIVITY	26
CROSS REACTIVITY	26
SAMPLE DILUTION	26
<b>14. REFERENCES</b>	<b>26</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>27</b>

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Tabelle 3: Wiederfindung n=2

Probe [ng/ml]	Spike [ng/ml]	slgA Erwartet [ng/ml]	slgA Gemessen [ng/ml]
103.7	150	253.7	279.7
103.7	75	178.7	194.7
103.7	50	153.7	158.7
103.7	25	128.7	141.2
100.3	150	250.3	272.4
100.3	75	175.3	212.9
100.3	50	150.3	165.4
100.3	25	125.3	126.5

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Lagerung/Storage: 2 - 8 °C

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of secretory IgA (slgA) in sputum and stool. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Secretory IgA (slgA) consists of two IgA monomers joined by the J-chain and furthermore a secretory component. It is secreted in plasma cells based in the lamina propria of mucosal membranes. Synthesis of slgA is independent from the synthesis of serum IgA. This means lack of serum IgA does not necessarily mean a lack of slgA. Secretory IgA is the major immunoglobulin in saliva, tears, colostrum, nasal mucous, mother's milk, tracheobronchial and gastrointestinal secretes. It plays a major role in preventing adherence of microorganisms to mucosal sites, in activation of the alternative complement pathway and in activating inflammatory reactions. Newborns are provided with slgA by mother's milk and are passively immunized against gastrointestinal infections.

### Indication

- Activity of the mucosa immune system (MIS)
- Proof of an imbalanced immunological barrier on the intestinal mucosa
- Autoimmune disease

### 3. PRINCIPLE OF THE TEST

In a first incubation step, the slgA in the samples is bound to polyclonal rabbit antibodies (in excess), which are immobilized to the surface of the microtitre wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labelled anti-slGA (POD-Antibody) antibody (polyclonal) is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of secretory IgA in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. Secretory IgA, present in the patient samples, is determined directly from this curve.

### 4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 8870MTP	one holder with precoated strips	12 x 8
K 8870WP	ELISA wash buffer concentrate 10x	1 x 100 ml
K 8870K	POD antibody, (polyclonal anti-slGA, Peroxidase-labeled), ready-to-use	1 x 15 ml
K 8870ST	Calibrators, lyophilized (0; 22.2; 66.6; 200; 600 ng/ml)	5 x 1 ml
K 8870KO	Control, lyophilized	1 x 1ml
K 8870TMB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 8870AC	ELISA stop solution, ready-to-use	1 x 7 ml

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized or bidistilled water
- Precision pipettes calibrated to deliver 10-1000  $\mu$ l
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Horizontal mixer
- Vortex-Mixer
- Microplate reader 450 nm

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** should be diluted with aqua dest. **1:10** before use (add 900 ml aqua bidest. to 100 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **calibrators and control** have to be reconstituted with **1000  $\mu$ l aqua bidest.** Allow the vial to stand for **10 minutes** and then mix thoroughly by gentle inversion to ensure complete reconstitution. Reconstituted calibrator and control are stable at 2-8 °C until the expiry date stated on the label.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2 – 8 °C.**

## 7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HCV or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The stop solution consists of diluted Sulfuric Acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date stated on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Saliva

To avoid variation we recommend to collect the sputum sample always at the same time of the day.

30 min before collection no food or liquid should be consumed. The samples are centrifuged at 3000 rpm for 10 min. Take the supernatant and dilute **1:2000** in ELISA wash buffer (add 10 µl in 1 ml and dilute this solution once more 50 µl in 1 ml - **100 µl** of this working solution is used).

### Faeces

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the ELISA wash buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. 1 ml supernatant is given into an eppendorf tube and centrifuged once more at 13.000 rpm for 5 min. The resulting supernatant could be stored at -20°C for about 1 month.

After thawing, the supernatant is centrifuged at 13.000 rpm for 2 min. Afterwards the supernatant is diluted **1:250** in ELISA wash buffer (4µl supernatant is added to 996 µl ELISA wash buffer). **100 µl** of this solution is used for the assay.

Immundiagnostik recommends the use of the sample tubes from Roche Diagnostics (Cat.-No. 745804) for sample preparation.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- Carry out the assay with the actual manual delivered with the kit.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** calibrators, control and patient samples (faeces and sputum diluted, see above).
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** Peroxidase-labelled anti-slgA antibody.
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** wash buffer.
7. Add **100 µl** TMB substrate solution.
8. Incubate for **5-10 minutes** at room temperature.


9. Add **50 µl** stop solution and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## 10. RESULTS

A calibration curve is constructed from the calibrators. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

Concentration [ng/ml]	600	200	66.6	22.2	0
OD mean value	2.209	1.286	0.692	0.310	0.118

These data are for demonstration only and cannot be used instead of data obtained from the actual assay.

## Saliva

For the calculation of the saliva values the results from the microplate reader has to be multiplied with **2.000**.

## Faeces

For the sIgA concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 250

dilution factor: 62,5 x 250 = 15625

Multiply the result with **15625** to get the real concentration. The dilution factor depends on the weight of the faeces.

## 11. LIMITATIONS

Samples with sIgA levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### *Expected values*

Secretory IgA (saliva) 102-471 µg/ml  
(Age >16 years)

Secretory IgA (faeces): 510-2040 µg/ml

### 13. PERFORMANCE CHARACTERISTICS

#### *Precision and reproducibility*

The precision (intra-assay variation) of the Immundiagnostik anti-slgA ELISA test was calculated from 20 replicate determinations on each of one samples.

Intra-Assay CV n= 20

Sample	slgA [ng/ml]	Intra-Assay CV [%]
1	77.7	5
2	92.5	9

The total precision (inter-assay variation) of the Immundiagnostik anti-slgA ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

Sample	slgA [ng/ml]	Inter-Assay CV [%]
1	102.4	8
2	1277.4	7.4

### Recovery

Two samples were spiked with slgA calibrator and measured with this assay.

Recovery n=2

Sample [ng/ml]	Spike [ng/ml]	slgA Expected [ng/ml]	slgA Measured [ng/ml]
103.7	150	253.7	279.7
103.7	75	178.7	194.7
103.7	50	153.7	158.7
103.7	25	128.7	141.2
100.3	150	250.3	272.4
100.3	75	175.3	212.9
100.3	50	150.3	165.4
100.3	25	125.3	126.5

*Sensitivity*

n=20

Sample	slgA Mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.171	0.099	13.4

*Cross reactivity*

No cross reactivity to other proteins in stool and saliva.

*Sample dilution*

Linearity n= 2

Two patient samples were diluted with ELISA wash buffer. The results are shown below:

Sample	Dilution	Expected [ng/ml]	Measured [ng/ml]
A	undiluted	126.8	126.8
	1:2	63.4	65.5
	1:4	31.7	35.1
	1:8	15.9	25.6
B	undiluted	184.9	184.9
	1:2	92.5	93.7
	1:4	46.2	52.1
	1:8	23.1	21.9

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## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from this.

10.10.2005 slgA\_10.10.2005.DOC